

# CT-520 Combined Table

- The new CT-520 Combined Table combines the modern design with the unique ability to effortlessly control the entire exam room environment.



## Features:

- Hydraulic pressure up and down movement with lower noise;
- Able to be equipped with arm, projector and so on;
- Economical background illumination, comfortable;
- Reasonable construction and clipper-built shape

## Specifications:

- Arm tilting angle:  $\pm 30^\circ$ ;
- Table rotate angle:  $90^\circ$ ;
- Table slip scope:  $\pm 150\text{mm}$ ;
- Chair rotate:  $\pm 60^\circ$ ;
- Chair height: 500-700mm;
- Weight: 200Kgs;
- Power: 120W, 220V, 50Hz/110V, 60Hz;
- Illumination bulb: 15W;
- Table size: 840mm $\times$ 405mm