

PM-300 Pachymeter



Features:

- High Light Large colour With touch screen
- Small, lightweight masterpiece
- User Friendly Interface
- 5 Programmable User Profiles
- Obtain fast and accurate Scans
- Echoes View
- USB and Mouse port
- Mass Data storage
- Built-in Printer
- Ultra-Portable

Specifications:

Corneal Thickness Measurement

Probe: 20MHz solid probe (non-gel type)

Tip Diameter: 1.5mm

Accuracy: 0.5µm

Measuring Range: 200-1300µm

Minimum Indicated Unit: 1µm

Measured Part: Corneal thickness up to 33 points can be memorized

Display

Large TFT colour Large touch screen

Memory

Stores thousands of data for eye

Printing Function

Built-in Graphic printer (Thermal type)

USB connection to U- Disk

USB and Mouse port

Power Requirements: AC100V~240V, 50 / 60Hz, 45VA Max

Dimensions / Weight: 240(W)×130(D)×205(H)mm / 1.7kg

